



Room Care R1

Toilet cleaner

Description

Ready-to-use, daily toilet bowl and urinal cleaner

Key properties

- Citric acid-based formulation
- Ideal viscosity for excellent adhesion
- Light, delicate and fresh fragrance
- Unique, patented O.N.T. air neutralising technology alters malodour molecules

Benefits

- Appropriate for daily use
- Effectively removes stains and deposits from toilet bowls and urinals
- Extended contact time on vertical/sloped surfaces for optimum cleaning efficiency
- Angled-neck squeezable bottle for fast and easy application
- Neutralises malodour molecules, enhancing the effect of the perfume
- Colour and alphanumeric coded to prevent application mistakes


Use instructions

Dosage:

Ready to use, no dilution required.

Application:

1. Flush the toilet/urinal and apply the product under the rim, around the bowl, onto water/flush paths and water outlets
2. Allow the product to act for 5 minutes
3. Scrub with a brush and flush to rinse
4. Increase the acting time (1 - 2 hours) in case of heavy calcium/lime build-up



© A.I.S.E.



Room Care R1

Technical data

Appearance: Clear, blue liquid

pH value (neat): < 2

pH value (in use): Relative density (20°C): ≈ 1.02

The above data is typical of normal production and should not be taken as a specification.

Safe handling and storage information

Full guidance on the handling and disposal of this product is provided in a separate Safety Data Sheet; sds.diverseym.com. Store in original closed containers away from extremes of temperature. Only for professional users / specialists.

Product compatibility

Do not use on acid-sensitive surfaces, e.g. enamel, marble, grouting. Immediately remove any drops or splashes from acid-sensitive surfaces and rinse with water. Do not mix with bleach or other cleaning products. Test material compatibility in a small, inconspicuous place before use.

Environmental information

The surfactants used in this product are biodegradable in compliance with the requirement of EU Detergent Regulation, EC 648/2004.

Available pack sizes

Room Care R1 is available in a 0.75L bottle, suitable for recycling.