



Room Care R2-plus

Hard surface cleaner and disinfectant - concentrate

Introduction to Room Care

Room Care is the unique globally available all encompassing solution for your guest room maintenance including chemicals, tools and methods to cover any room cleaning needs, adapted to any hotel size and housekeeping structure.

Superior formulations adapted to deal effectively and safely with all the cleaning and disinfections requirements of housekeeping, easy to recognise product names as well as colour coding throughout the system, a set of dedicated tools, accessories and machines make Room Care the perfect choice.

All Room Care products come with the same 'Fresh Linen' signature fragrance designed for the lodging industry. It utilises the O.N.T. technology for effective malodour neutralisation and long lasting fragrance.

Diversey's expertise, environmental commitment, and certification against ISO 9001, ISO 14001 and the AISE Charter for Sustainable Cleaning provide further guarantee of quality and efficacy.

Product description

Room Care R2-plus is a concentrated cleaner and disinfectant for the hygienic cleaning of all water resistant hard surfaces. To be used only with Divermite/DQFM dosing and dilution system.

Key properties

- Formulated for the hygienic cleaning of telephones and other hard, water-resistant surfaces
- Slightly alkaline
- Non-abrasive
- Light, delicate and fresh fragrance
- Unique, patented ONT air neutralising technology alters malodour molecules

Benefits

- Leaves surfaces hygienically clean
- Quickly and gently removes oily and water-soluble residues
- Regular use hinders lime scale build-up
- Leaves a long-lasting pleasant scent behind
- Neutralises malodour molecules, enhancing the effect of the perfume
- Colour and alphanumeric coded to prevent application mistakes





Room Care R2-plus

Use instructions

When using with Divermite Room Care Dosing Unit, fill keyed Room Care R2-plus 300ml application bottle with water to the fill line. Push the Divermite dispenser twice to dispense 30mls of product. Securely screw on the trigger spray. Ensure the nozzle is closed and gently shake the bottle to mix the product.

When using with a DQFM dilution dispenser unit, fill the Room Care R2-plus keyed spray bottle (300ml or 750ml) to the fill level mark directly from the DQFM dispenser. Securely screw-on the spray trigger.

Application

Spray a small quantity onto a damp cloth/sponge pad and clean the surface. Rinse with fresh water or wipe with a clean, damp cloth. Buff with a dry cloth for extra shine.

NB: Between refills, rinse out bottle with clean water.

Product Compatibility

Test material compatibility in a small, inconspicuous place before use. Do not use with other dilution systems or methods. Do not use undiluted.

Technical data

Appearance	Clear, green liquid	pH-value neat	10.5 - 11.0
Relative density (20°C)	1.030	pH-value in use	10.5 +/- 0.5 at 10% solution

The above data is typical of normal production and should not be taken as a specification.

Product	Pack size	Article code
Room Care R2-plus – Hygienic hard surface cleaner	2 x 1.5 L pouch	100941228
Bottle Kit (300ml)	6x0.3L	1204321
Bottle Kit (750ml)	6x0.75L	1204315

Safe handling and storage information

Rinse and dry hands after use. The use of a skin reconditioning cream, such as Soft Care™ Derasoft, may be beneficial. For prolonged contact, protection for the skin may be necessary. Do not breathe spray. Use only in well-ventilated areas. Full guidance on the handling and disposal of this product is provided in a separate Safety Data Sheet. Only for professional users/specialists. Store in original closed containers away from extremes of temperature.

Environmental information

The surfactants used in this product are biodegradable in compliance with the requirement of EC Directives 73/404/EEC and 73/405/EEC and their subsequent amendments.

Microbiological data

Room Care R2-plus passed EN 1276 at 10% dilution in water of standardised hardness and in dirty conditions (0.3% bovine albumin) at a contact time of 1 minute.

Room Care R2-plus passed EN 13697 at 10% dilution in water of standardised hardness and in dirty conditions (0.3% bovine albumin) for referenced strains *Staphylococcus aureus*, *Enterococcus hirae*, *Escherichia coli* and *Pseudomonas aeruginosa* at a contact time of 5 minutes. Yeasticidal activity in 15 minutes when diluted at 10% for the referenced strain *Candida albicans*.