

PRODUCT SPECIFICATION

Product name	Laundry Bag
Range	Premier
Ref	2845
Marking	Product name, Premier Logo, Code Ref
Material	Polythene with degradable plastic additive, starch based soluble seam and tie closure
Colour	Red
Sterilisation	Non-Sterile



Primary use	Intended to provide a hygienic way of isolating and transporting dirty or soiled linen.
--------------------	---

Size (Open)	Size (Closed)	Thickness	Packaging
711x762mm (28x30")	460x762mm (18x30")	20mu	200 Bags per carton

Storage conditions	Store in a cool dry place away from heat and direct sunlight
Additional information	<p>The product is latex free and non-sterile.</p> <p>Load the bag directly into the washing machine to eliminate the need to open bags to sort and/or load the contaminated linen as well as the potential exposure to the handlers.</p> <p>During the agitation cycle, the water-soluble seam and tie closure dissolve, releasing the contents of the bag into the washing solution.</p> <p>At the conclusion of the wash cycle, the used and disinfected portion of the polyethylene bag is removed and discarded or recycled.</p> <p>Minimum operating temperature for soluble strip is 23°C.</p> <p>This product is not suitable for tunnel washers and tumble-drying machines.</p>

REVISION:	4.0	DATE	05/01/2021	AUTHOR	A Gidley
------------------	-----	-------------	------------	---------------	----------

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare T: +44 (0) 1530 278 111 W: www.365healthcare.com	Unit 6, Delta Park Industrial Estate Millmarsh Lane, London Enfield, EN3 7QJ
---	--